

Product Datasheet

IL-10 antibody (orb433206)

Description	Mouse monoclonal antibody to IL-10
Species/Host	Mouse
Reactivity	Bovine
Conjugation	Unconjugated
Tested Applications	ELISA, ELISpot, FC
Immunogen	Plasmid cDNA encoding bovine IL-10.
Target	IL-10
Preservatives	0.09% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.