

Product Datasheet

Granzyme B antibody (orb433218)

Description	Mouse monoclonal antibody to Granzyme B
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant Granzyme B.
Target	GRANZYME B
Preservatives	0.09% Sodium Azide (NaN ₃). TRIS buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.