

Product Datasheet

BETA AMYLOID antibody (orb433270)

Description	Mouse monoclonal antibody to BETA AMYLOID
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, IP
Immunogen	Synthetic peptide consisting of residues 1 - 17 of the amyloid beta protein.
Target	BETA AMYLOID
Preservatives	0.09% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	(aa1-17)
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.