

Product Datasheet

CD195 antibody (FITC) (orb433274)

Description	Rat monoclonal antibody against CD195 conjugated to FITC
Species/Host	Rat
Reactivity	Human
Conjugation	FITC
Tested Applications	FC
Immunogen	CHO cells transfected with human CCR5
Target	CD195
Preservatives	0.09% sodium azide (NaN ₃)1% bovine serum albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	0.1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.