

Product Datasheet

MHC CLASS I antibody (orb433312)

Description	Mouse monoclonal antibody to MHC CLASS I
Species/Host	Mouse
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC, IP
Immunogen	Purified H-2Kb and H-2Db MHC-I molecules.
Target	MHC CLASS I
Preservatives	0.09% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.