

Product Datasheet

MHC CLASS I antibody (FITC) (orb433313)

Description	Mouse monoclonal antibody against MHC CLASS I conjugated to FITC
Species/Host	Mouse
Reactivity	Mouse
Conjugation	FITC
Tested Applications	FC
Immunogen	Purified H-2Kb and H-2Db MHC-I molecules.
Target	MHC CLASS I
Preservatives	0.09% sodium azide (NaN ₃)1% bovine serum albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	0.1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.