

Product Datasheet

5-METHYLCYTIDINE antibody (orb433324)

Description	Mouse monoclonal antibody to 5-METHYLCYTIDINE
Species/Host	Mouse
Reactivity	Aves
Conjugation	Unconjugated
Tested Applications	DOT, ELISA, FACS, ICC, IF, IHC-Fr, IHC-P, IP
Target	5-METHYLCYTIDINE
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.