

Product Datasheet

CD79b antibody (orb433327)

Description	Mouse monoclonal antibody to CD79b
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IP, WB
Immunogen	Peptide corresponding to an extracellular region of CD79 beta
Target	CD79b
Preservatives	0.09% sodium azide (NaN ₃). TRIS buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.