

Product Datasheet

CD45 antibody (Biotin) (orb433345)

Description	Mouse monoclonal antibody against CD45 conjugated to Biotin
Species/Host	Mouse
Reactivity	Sheep
Conjugation	Biotin
Tested Applications	FC
Immunogen	Ovine efferent lymphatic duct lymphocytes.
Target	CD45
Preservatives	0.09% sodium azide (NaN ₃)1% bovine serum albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to biotin - liquid
Concentration	0.1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.