

Product Datasheet

CD45 antibody (RPE) (orb433348)

Description	Mouse monoclonal antibody against CD45 conjugated to RPE
Species/Host	Mouse
Reactivity	Sheep
Conjugation	RPE
Tested Applications	FC
Immunogen	Ovine efferent lymphatic duct lymphocytes.
Target	CD45
Preservatives	0.09% sodium azide (NaN ₃)1% bovine serum albumin5% sucrose. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.