

Product Datasheet

MHC CLASS II DQ DR POLYMORPHIC antibody (RPE) (orb433353)

Description	Mouse monoclonal antibody against MHC CLASS II DQ DR POLYMORPHIC conjugated to RPE
Species/Host	Mouse
Reactivity	Sheep
Conjugation	RPE
Tested Applications	FC
Immunogen	Ovine alveolar macrophages.
Target	MHC CLASS II DQ DR POLYMORPHIC
Preservatives	0.09% sodium azide (NaN ₃)1% bovine serum albumin5% sucrose. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.