

Product Datasheet

MHC CLASS II DQ DR POLYMORPHIC antibody (FITC) (orb433356)

Description	Mouse monoclonal antibody against MHC CLASS II DQ DR POLYMORPHIC conjugated to FITC
Species/Host	Mouse
Reactivity	Sheep
Conjugation	FITC
Tested Applications	FC
Target	MHC CLASS II DQ DR POLYMORPHIC
Preservatives	0.09% sodium azide (NaN ₃)1% bovine serum albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	0.1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.