

Product Datasheet

CD206 antibody (orb433360)

Description	Rat monoclonal antibody to CD206
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FACS, IHC-Fr, IP
Immunogen	Chimaeric CRD4-7-Fc protein
Target	CD206
Preservatives	None Present. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Low Endotoxin
Expiration Date	12 months from date of receipt.