

Product Datasheet

CD160 antibody (orb433374)

Description	Mouse monoclonal antibody to CD160
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IP
Immunogen	Human NK cell line YT2C2.
Target	CD160
Preservatives	0.09% sodium azide (NaN ₃)1% bovine serum albumin. Phosphate buffered saline
Form/Appearance	Purified IgM - liquid
Concentration	0.1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgM
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.