

Product Datasheet

LYMPHOTOXIN BETA RECEPTOR antibody (RPE) (orb433393)

Description	Rat monoclonal antibody against LYMPHOTOXIN BETA RECEPTOR conjugated to RPE
Species/Host	Rat
Reactivity	Mouse
Conjugation	RPE
Tested Applications	FC
Immunogen	Recombinant extracellular domain of LT beta R.
Target	LYMPHOTOXIN BETA RECEPTOR
Preservatives	0.09% sodium azide (NaN ₃)1% bovine serum albumin5% sucrose. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.