

Product Datasheet

RITUXIMAB antibody (orb433424)

Description	Rat monoclonal antibody to RITUXIMAB
Species/Host	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Immunogen	F(ab)2 fragment of Rituximab
Target	RITUXIMAB
Preservatives	0.1% sodium azide (NaN ₃). TRIS buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.