

Product Datasheet

SLA CLASS I antibody (orb433428)

Description	Mouse monoclonal antibody to SLA CLASS I
Species/Host	Mouse
Reactivity	Porcine
Conjugation	Unconjugated
Tested Applications	FACS, IP
Immunogen	Porcine peripheral blood mononuclear cells.
Target	SLA CLASS I
Preservatives	0.09% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.