

Product Datasheet

ESTROGEN RECEPTOR BETA 2 antibody (orb433442)

Description	Mouse monoclonal antibody to ESTROGEN RECEPTOR BETA 2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Synthetic peptide derived from the C terminus of the human estrogen receptor beta 2 isomer, CGMKMETLLPEATMEQ.
Target	ESTROGEN RECEPTOR BETA 2
Preservatives	0.1% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.