

Product Datasheet

CD147 antibody (orb433460)

Description	Rat monoclonal antibody to CD147
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC, IH-Fr, WB
Immunogen	Fusion protein mCD147dO-2rCD4d3 +4
Target	CD147
Preservatives	0.09% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.