

## Product Datasheet

### CD147 antibody (orb433461)

|                            |                                                                   |
|----------------------------|-------------------------------------------------------------------|
| <b>Description</b>         | Rat monoclonal antibody to CD147                                  |
| <b>Species/Host</b>        | Rat                                                               |
| <b>Reactivity</b>          | Mouse                                                             |
| <b>Conjugation</b>         | Unconjugated                                                      |
| <b>Tested Applications</b> | FC, IH-Fr, WB                                                     |
| <b>Immunogen</b>           | Fusion protein mCD147dO-2rCD4d3 +4                                |
| <b>Target</b>              | CD147                                                             |
| <b>Preservatives</b>       | 0.09% sodium azide (NaN <sub>3</sub> ). Phosphate buffered saline |
| <b>Form/Appearance</b>     | Purified IgG - liquid                                             |
| <b>Concentration</b>       | 1.0 mg/ml                                                         |
| <b>Storage</b>             | +4°C, -20°C if preferred                                          |
| <b>Note</b>                | For research use only                                             |
| <b>Isotype</b>             | IgG1                                                              |
| <b>Clonality</b>           | Monoclonal                                                        |
| <b>Purity</b>              | Purified                                                          |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                   |