

Product Datasheet

DECTIN-1 antibody (orb433471)

Description	Rat monoclonal antibody to DECTIN-1
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC, IH-Fr, IP
Immunogen	Dectin-1 transfected NIH3T3 cells and recombinant soluble Dectin-1.
Target	DECTIN-1
Preservatives	0.09% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.