

Product Datasheet

DECTIN-1 antibody (FITC) (orb433477)

Description	Rat monoclonal antibody against DECTIN-1 conjugated to FITC
Species/Host	Rat
Reactivity	Mouse
Conjugation	FITC
Tested Applications	FC
Immunogen	Dectin-1 transfected NIH3T3 cells and recombinant soluble Dectin-1.
Target	DECTIN-1
Preservatives	0.09% sodium azide (NaN ₃)1% bovine serum albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	0.5 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.