

Product Datasheet

CD107b antibody (orb433488)

Description	Rat monoclonal antibody to CD107b
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC, IH-Fr, IHC-P, IP, WB
Immunogen	Glycoproteins purified from mouse peritoneal macrophage membranes.
Target	CD107b
Preservatives	0.09% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.