

Product Datasheet

CD29 antibody (orb433501)

Description	Hamster monoclonal antibody to CD29
Species/Host	Hamster
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FA, FACS, IP
Immunogen	Purified mouse VLA-4 antigen.
Target	CD29
Preservatives	None present. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	Low Endotoxin
Expiration Date	12 months from date of receipt.