

Product Datasheet

CD29 antibody (RPE) (orb433506)

Description	Hamster monoclonal antibody against CD29 conjugated to RPE
Species/Host	Hamster
Reactivity	Mouse
Conjugation	RPE
Tested Applications	FC
Immunogen	Purified mouse VLA-4 antigen.
Target	CD29
Preservatives	0.09% sodium azide (NaN ₃)1% bovine serum albumin5% sucrose. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.