

Product Datasheet

CD16/CD32 antibody (orb433512)

Description	Rat monoclonal antibody to CD16/CD32
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	PU5 1.8 IOE7 Balb/c mouse cell line.
Target	CD16/CD32
Preservatives	0.09% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.