

Product Datasheet

CD312 antibody (orb433554)

Description	Mouse monoclonal antibody to CD312
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IH-Fr, IP, WB
Immunogen	NIH-3T3 cells stably transfected with EMR2 (EGF1-5) cDNA.
Target	CD312
Preservatives	0.09% Sodium Azide. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.