

Product Datasheet

DAZL antibody (orb433564)

Description	Mouse monoclonal antibody to DAZL
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Immunogen	Synthetic peptide corresponding to sequence within the C terminal domain of human DAZL (CRVHHFRRSRAMLKSV).
Target	DAZL
Preservatives	0.1% Sodium Azide (NaN ₃)
Form/Appearance	Tissue Culture Supernatant - liquid
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	S/N
Expiration Date	12 months from date of receipt.