

Product Datasheet

CD223 antibody (orb433612)

Description	Rat monoclonal antibody to CD223
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FA, FACS
Immunogen	Murine CD223 Ig fusion protein.
Target	CD223
Preservatives	None present. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Low Endotoxin
Expiration Date	12 months from date of receipt.