

Product Datasheet

Gr-1 antibody (RPE) (orb433620)

Description	Rat monoclonal antibody against Gr-1 conjugated to RPE
Species/Host	Rat
Reactivity	Mouse
Conjugation	RPE
Tested Applications	FC
Immunogen	Normal murine bone marrow cells.
Target	Gr-1
Preservatives	0.09% sodium azide (NaN ₃)1% bovine serum albumin5% sucrose. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.