

Product Datasheet

Gr-1 antibody (orb433622)

Description	Rat monoclonal antibody to Gr-1
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FACS, IHC-Fr, IP, WB
Immunogen	Normal murine bone marrow cells.
Target	Gr-1
Preservatives	0.09% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.