

Product Datasheet

ER-TR7 antibody (orb433646)

Description	Rat monoclonal antibody to ER-TR7
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FACS, ICC, IF, IHC-Fr
Immunogen	Isolated C3H thymic stromal cells.
Target	ER-TR7
Preservatives	0.09% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.