

Product Datasheet

CD34 antibody (Biotin) (orb433665)

Description	Mouse monoclonal antibody against CD34 conjugated to Biotin
Species/Host	Mouse
Reactivity	Canine
Conjugation	Biotin
Tested Applications	FC
Immunogen	Canine CD34 fusion protein.
Target	CD34
Preservatives	0.09% Sodium Azide 1% Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Biotin - liquid
Concentration	0.1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.