

Product Datasheet

IFN ALPHA antibody (orb433669)

Description	Mouse monoclonal antibody to IFN ALPHA
Species/Host	Mouse
Reactivity	Gallus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant chicken Interferon alpha.
Target	IFN ALPHA
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	0.5 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.