

Product Datasheet

CD45 antibody (RPE) (orb433671)

Description	Mouse monoclonal antibody against CD45 conjugated to RPE
Species/Host	Mouse
Reactivity	Gallus
Conjugation	RPE
Tested Applications	FC
Immunogen	COS cells transfected with CD45 gene.
Target	CD45
Preservatives	0.09% Sodium Azide 1% Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.