

Product Datasheet

IgA antibody (HRP) (orb433705)

Description	Mouse monoclonal antibody against IgA conjugated to HRP
Species/Host	Mouse
Reactivity	Bovine
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Purified bovine Immunoglobulin
Target	IgA
Preservatives	0.01%Thiomersal. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.