

Product Datasheet

IgM antibody (HRP) (orb433714)

Description	Mouse monoclonal antibody against IgM conjugated to HRP
Species/Host	Mouse
Reactivity	Bovine
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Bovine peripheral blood mononuclear cells depleted of T lymphocytes.
Target	IgM
Preservatives	HRP Stabiliser (BUF052A)0.01% Thiomersal. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Concentration	0.5 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.