

## Product Datasheet

### CD88 antibody (Biotin) (orb433727)

|                            |                                                                       |
|----------------------------|-----------------------------------------------------------------------|
| <b>Description</b>         | Rat monoclonal antibody against CD88 conjugated to Biotin             |
| <b>Species/Host</b>        | Rat                                                                   |
| <b>Reactivity</b>          | Mouse                                                                 |
| <b>Conjugation</b>         | Biotin                                                                |
| <b>Tested Applications</b> | FC                                                                    |
| <b>Immunogen</b>           | RBL-2H3 transfected cells expressing murine C5aR.                     |
| <b>Target</b>              | CD88                                                                  |
| <b>Preservatives</b>       | 0.09% Sodium Azide 1% Bovine Serum Albumin. Phosphate buffered saline |
| <b>Form/Appearance</b>     | Purified IgG conjugated to Biotin - liquid                            |
| <b>Concentration</b>       | 0.1 mg/ml                                                             |
| <b>Storage</b>             | +4°C, -20°C if preferred                                              |
| <b>Note</b>                | For research use only                                                 |
| <b>Isotype</b>             | IgG2a                                                                 |
| <b>Clonality</b>           | Monoclonal                                                            |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                       |