

## Product Datasheet

### CD88 antibody (orb433730)

|                            |                                                   |
|----------------------------|---------------------------------------------------|
| <b>Description</b>         | Rat monoclonal antibody to CD88                   |
| <b>Species/Host</b>        | Rat                                               |
| <b>Reactivity</b>          | Mouse                                             |
| <b>Conjugation</b>         | Unconjugated                                      |
| <b>Tested Applications</b> | FC                                                |
| <b>Immunogen</b>           | RBL-2H3 transfected cells expressing murine C5aR. |
| <b>Target</b>              | CD88                                              |
| <b>Preservatives</b>       | 0.09% Sodium Azide. Phosphate buffered saline     |
| <b>Form/Appearance</b>     | Purified IgG - liquid                             |
| <b>Concentration</b>       | 1.0 mg/ml                                         |
| <b>Storage</b>             | +4°C, -20°C if preferred                          |
| <b>Note</b>                | For research use only                             |
| <b>Isotype</b>             | IgG2b                                             |
| <b>Clonality</b>           | Monoclonal                                        |
| <b>Purity</b>              | Purified                                          |
| <b>Expiration Date</b>     | 12 months from date of receipt.                   |