

Product Datasheet

CD138 antibody (orb433741)

Description	Mouse monoclonal antibody to CD138
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FACS, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	U266 cell line.
Target	CD138
Preservatives	0.09% Sodium Azide 1% Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	0.1mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.