

Product Datasheet

CD279 antibody (orb433758)

Description	Rat monoclonal antibody to CD279
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	PD-1 transferred BHK cells.
Target	CD279
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.