

Product Datasheet

GLUCOCORTICOID RECEPTOR antibody (orb433766)

Description	Mouse monoclonal antibody to GLUCOCORTICOID RECEPTOR
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	26 amino acid peptide corresponding to residues 150-176 on human GCR linked to thyroglobulin.
Target	GLUCOCORTICOID RECEPTOR
Preservatives	0.1% Sodium Azide (NaN ₃) 0.1% Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.