

Product Datasheet

PHOSPHOTYROSINE antibody (Biotin) (orb433770)

Description	Mouse monoclonal antibody against PHOSPHOTYROSINE conjugated to Biotin
Species/Host	Mouse
Reactivity	Chemical
Conjugation	Biotin
Tested Applications	ELISA, ICC, IF, WB
Target	PHOSPHOTYROSINE
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Biotin - liquid
Concentration	0.5 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.