

Product Datasheet

PHOSPHOTYROSINE antibody (HRP) (orb433771)

Description	Mouse monoclonal antibody against PHOSPHOTYROSINE conjugated to HRP
Species/Host	Mouse
Reactivity	Chemical
Conjugation	HRP
Tested Applications	ELISA
Target	PHOSPHOTYROSINE
Preservatives	0.01% Thiomersal50% Glycerol. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.