

Product Datasheet

CD71 antibody (orb433776)

Description	Rat monoclonal antibody to CD71
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FA, FACS, IHC-Fr
Immunogen	Mouse transformed endothelioma cell line t.end1.
Target	CD71
Preservatives	None present. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Low Endotoxin
Expiration Date	12 months from date of receipt.