

Product Datasheet

IgY antibody (orb433785)

Description	Mouse monoclonal antibody to IgY
Species/Host	Mouse
Reactivity	Aves
Conjugation	Unconjugated
Tested Applications	ELISA, FC, WB
Immunogen	Purified Pekin duck yolk IgY.
Target	IgY
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	LIGHT CHAIN
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.