


Product Datasheet

Coiled-coil domain-containing protein 106 Antibody (orb43380)

Description	Coiled-coil domain-containing protein 106 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q9BWC9
Expiration Date	12 months from date of receipt.



Western blot analysis of Coiled-coil dom...