

Product Datasheet

Zinc finger CCHC domain-containing protein 10 Antibody (orb43381)

Description Zinc finger CCHC domain-containing protein 10
Antibody

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications IF, IHC, WB


Storage Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Note For research use only

Clonality Recombinant

Uniprot ID [Q5EB97](https://www.uniprot.org/entry/Q5EB97)

Expiration Date 12 months from date of receipt.



Western
blot
analysis of
Zinc finger
CCH...