

Product Datasheet

CD74 antibody (RPE) (orb433815)

Description	Mouse monoclonal antibody against CD74 conjugated to RPE
Species/Host	Mouse
Reactivity	Human
Conjugation	RPE
Tested Applications	FC
Immunogen	B lymphoblastoid cell line HFB1.
Target	CD74
Preservatives	0.09% Sodium Azide (NaN ₃)1% Bovine Serum Albumin5% Sucrose. Phosphate buffered saline
Form/Appearance	Purified Protein conjugated to R. Phycoerythrin (RPE) - lyophilized
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.