

Product Datasheet

CD172a antibody (orb433836)

Description	Mouse monoclonal antibody to CD172a
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FACS, IHC-Fr
Immunogen	Monocyte-derived dendritic cells.
Target	CD172a
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.